

Comparison of six different calprotectin assays for the assessment of inflammatory bowel disease

Delphine Labaere¹, Annick Smismans¹, August Van Olmen², Paul Christiaens², Geert D'Haens^{2,3}, Veerle Moons², Pieter-Jan Cuyle², Johan Frans¹ and Peter Bossuyt²

Abstract

Background and objectives: Faecal calprotectin is a valuable noninvasive marker for inflammatory bowel disease (IBD). The aim of our study was to determine the correlation between six different calprotectin assays and compare their performance for diagnosis and follow up of IBD.

Methods: Thirty-one patients with suspected IBD and 31 patients in follow up were included. We determined calprotectin by means of three rapid immunochromatographic tests, two enzyme-linked immunosorbent assays, and one automated fluoroimmunoassay. Results were correlated with endoscopic and histological findings.

Results: Although all methods correlated significantly, slopes and intercepts differed extensively, with up to 5-fold quantitative differences between assays. Sensitivity and specificity for diagnosis of IBD were 82–83 and 84–89%, respectively. For follow up, sensitivity in detecting mild ulcerative colitis was 71–100%. In moderate-to-severe ulcerative colitis, sensitivity was 100% for all assays. Specificity was 67–86% in both subgroups. In Crohn's disease, only moderate-to-severe disease could be differentiated from remission, with sensitivity 83–86% and specificity 75% for all tests.

Conclusions: All calprotectin assays showed comparable clinical performance for diagnosis of IBD. For follow up, performance was acceptable, except for mild Crohn's disease. Because of the large quantitative differences, further efforts are needed to standardize calprotectin assays.

Keywords

Calprotectin, Crohn's disease, diagnosis, follow up, inflammatory bowel disease, ulcerative colitis

Received: 30 October 2013; accepted: 3 December 2013

Introduction

Ileocolonoscopy with histopathological sampling is generally considered the gold standard in the diagnosis and follow up of inflammatory bowel disease (IBD).¹ However, endoscopy procedures and the required bowel preparation are expensive, time consuming, and inconvenient to the patient.

Several biochemical parameters in plasma and stool have been investigated for their ability to detect bowel inflammation, of which faecal calprotectin (FC) has shown the best test performance until now.^{2,3} Calprotectin is a protein that is released by neutrophils in the bowel in case of inflammation, as is the case in IBD. In patients suspected for IBD, a negative FC indicates a low chance of having bowel inflammation (high negative predictive value), whereas an elevated FC

needs further investigation with endoscopy and histology. A recent meta-analysis calculated that calprotectin screening might reduce the need for endoscopic procedures in adults with suspected IBD in up to 67%.⁴

In patients with known IBD, FC is high when active inflammation is present and low in case of remission. Recent promising results were obtained with FC to

¹Department of Laboratory Medicine, Imelda General Hospital, Bonheiden, Belgium

²Department of Gastroenterology, Imelda General Hospital, Bonheiden, Belgium

³Department of Gastroenterology, Academic Medical Center, Amsterdam, The Netherlands

Corresponding author:

Delphine Labaere, Department of Laboratory Medicine, Imelda Hospital, Imeldalaan 9, 2820 Bonheiden, Belgium.
Email: delphine.labaere@gmail.com

monitor disease activity in IBD patients, performing better than traditional inflammatory markers such as C-reactive protein.^{2,5}

Calprotectin is often measured with enzyme-linked immunosorbent assays (ELISA). Since ELISA is time consuming and mostly suited for analysing samples in batch, faster and more user-friendly techniques have been developed. However, very few data are available on the comparability of these different methods. The goal of this study was to evaluate six different calprotectin tests and to determine their performance for diagnosis and follow up of IBD patients.

Materials and methods

Study population

Sixty-four adult patients (≥ 16 years old), referred for colonoscopy to the gastroenterology department of our hospital, provided a stool sample for calprotectin measurement. Patients were prospectively included based on the indication for the endoscopy, being either clinical suspicion of IBD or the follow up of known IBD. Demographic characteristics of the study population are shown in Table 1. The patients with suspected IBD ($n=33$) were selected based on one or more of the following criteria: chronic diarrhoea and/or abdominal pain, iron deficiency anaemia, unexplained weight loss, or a family history of IBD. Two of these patients had to be excluded: in one case, the histopathological and endoscopic evaluation was inconclusive; in the second patient, the consistency of the stool was extremely watery and calprotectin was below the detection limit of all assays. This resulted in a population of 31 included patients with suspected IBD. The study population of known IBD patients undergoing follow-up endoscopy consisted of 31 patients, of which 17 had ulcerative colitis (UC) and 14 had Crohn's disease (CD). There were no exclusion criteria concerning medical history or drug intake, because this information was not available for most of our study patients.

Colonoscopy

Bowel preparation was done at the outpatient clinic with electrolyte or polyethylene glycol solutions.

Table 1. Demographic characteristics of the study populations

Study group	n	M/F	Age (years)		
			Median	Range	IQR
Suspected IBD	31	10/21	36	16–75	30–45
Follow up IBD	31	17/14	50	24–83	42–60

F, female; IBD, inflammatory bowel disease; M, male.

Ileocolonoscopies were performed by five experienced gastroenterologists. Endoscopic disease activity was evaluated by means of the Mayo endoscopic subscore⁶ for UC and the simple endoscopic score for CD (SES-CD).⁷

During endoscopy, biopsy specimens targeted at the most severely affected areas and from the edge of ulcers were obtained. If no lesions were present, biopsies were collected from random sites. Routine histology was performed by three experienced pathologists.

Calprotectin measurement

From each patient, a faecal sample was obtained on the day of the endoscopy. To minimize interference by the bowel preparation solution, the first stool that was obtained after starting the bowel preparation was collected. Sample collection was coordinated by the nurses of the outpatient clinic.

We determined FC by means of three rapid immunochromatographic assays (Bühlmann Quantum Blue, Eurospital Calfast, Biotec Certest), two ELISAs (Eurospital and Calprolab), and one automated immunoassay (Phadia EliA). The tests' characteristics are shown in Table 2.

The rapid tests are based on lateral flow immunochromatography, resulting in a line on the test strip if the analyte is present in the sample. Bühlmann Quantum Blue and Eurospital Calfast produce quantitative results, as the density of the test line is measured by an automated reader. Biotec Certest is a semi-quantitative test, with visual interpretation of the test line(s), resulting in a value <50 , $50\text{--}200$, or >200 $\mu\text{g/g}$. ELISA is a quantitative method on a microtitre plate. In our study, ELISAs were performed manually and optical densities were read afterwards by a PhD plate reader (BioRad). EliA from Phadia is a recently developed quantitative enzyme fluoroimmunoassay performed on an automated analyser (Immunocap 250).

All analyses were performed in the laboratory of the Imelda Hospital strictly according to the manufacturer's instructions. Stool samples were stored at $2\text{--}8^\circ\text{C}$ for a maximum of 1 week. Afterwards, stool extraction was completed using commercial extraction devices recommended by the manufacturer, as shown in Table 2. Since each assay has its own extraction buffer and extraction procedure, it was not possible to assess the performance of all assays using a common faecal extract. Extracts were centrifuged and stored for a maximum of 3 months at -20°C before analysis. Only Biotec Certest had to be performed immediately after extraction, according to the manufacturer's instructions.

Three samples of the diagnostic test group could not be analysed with Eurospital Calfast and Calprest assays, because of the limited amount of available test kits.

Statistics

Correlation between the different assays was calculated using nonparametric Passing Bablok regression analysis. Only those results within the measurement range of the assay were included. The degree of correlation between assays was determined by the nonparametric Spearman rank correlation coefficient. The significance level was determined by the associated *p*-value set at <0.05.

Receiver operating characteristic (ROC) analyses were performed to determine the optimal cut offs that provided the highest sum of sensitivity and specificity. The following performance characteristics were calculated: sensitivity, specificity, positive and negative predictive values, and positive and negative likelihood ratios.

All statistical analyses and graphics were performed using GraphPad Prism version 4, except the Passing Bablok regression analysis which was performed using Medcalc Software version 12.

Ethical considerations

All patients participated at a voluntary basis and provided informed consent. The study was approved by the local ethics committee of the Imelda Hospital (registration number B689201112774).

Results

Correlation between assays

Method comparison was performed for the five quantitative assays: Bühlmann Quantum Blue, Calprolab, Eurospital Calfast, Eurospital Calprest, and

Phadia EliA. The Passing Bablok regression and Spearman rank analysis characteristics are summarized in Table 3. Significant correlations (*p*<0.05) between all assays were obtained. Correlation coefficients ranged from 0.65 to 0.93. However, the slopes and/or intercepts differed extensively, with up to 5-fold quantitative differences between assays.

Clinical performance

Diagnosis. For the 31 enrolled patients with suspected IBD, a diagnosis of UC and CD was confirmed by ileocolonoscopy and histology in five and seven

Table 3. Passing Bablok analysis and Spearman rank correlation between the different assays

Method X	Method Y	<i>n</i>	Passing Bablok regression equation	Spearman rank correlation coefficient	<i>p</i> -value
Calprest	QB BM	37	$y=2.39x-19$	0.91	<0.001
Calprest	Phadia EliA	24	$y=1.67x-80$	0.70	<0.001
Calprest	Calprolab	29	$y=2.91x-108$	0.85	<0.001
Calprest	Calfast	27	$y=0.76x+19$	0.82	<0.001
Calfast	QB BM	31	$y=5.37x-266$	0.69	<0.001
Calfast	Phadia EliA	26	$y=1.66x-71$	0.65	<0.001
Calfast	Calprolab	31	$y=4.85x-227$	0.75	<0.001
Calprolab	QB BM	33	$y=1.12x+51$	0.93	<0.001
Calprolab	Phadia EliA	32	$y=0.31x+6$	0.65	<0.001
QB BM	Phadia EliA	28	$y=0.34x-15$	0.77	<0.001

n, number of samples available for regression analysis (results within measurement range of both tests); QB BM, Quantum Blue Bühlmann

Table 2. Characteristics of the six tested calprotectin assays

Assay (manufacturer)	Method	Extraction device	Measurement range (µg/g)
Calprest (Eurospital, Trieste, Italy)	ELISA	Eurospital extraction device	15–500
Calprolab (Calpro, Oslo, Norway)	ELISA	Roche faecal extraction device	25–2500
Calprotectin EliA (Phadia/Thermo-Fisher, Uppsala, Sweden)	Enzyme fluoroimmunoassay	Roche faecal extraction device	15–3000
Calfast (Eurospital, Trieste, Italy)	Quantitative immunochromatography (automated reading)	Eurospital extraction device	50–300
Calprotectin Quantum Blue (Bühlmann, Schönenbuch, Switzerland)	Quantitative immunochromatography (automated reading)	Roche faecal extraction device	30–300
Certest (Biotec, Zaragoza, Spain)	Semi-quantitative immunochromatography (visual reading)	Certest extraction device	50–200

ELISA, enzyme-linked immunosorbent assay.

Table 4. Diagnostic accuracy of the different calprotectin assays

	Cut off ($\mu\text{g/g}$)	IBD (n)	Non-IBD (n)	Median FC in non-IBD ($\mu\text{g/g}$)	Median FC in IBD ($\mu\text{g/g}$)	Ratio median FC IBD/non-IBD ^a	Sensitivity	Specificity	PPV	NPV	LR+	LR-
Calprolab	50	12	19	<25	239	10	83	89	83	89	7.9	0.19
QB BM	50	12	19	<30	425	14	83	68	63	87	2.7	0.24
	75						83	84	77	89	5.28	0.20
Calfast	70	11	17	<50	132	3	82	88	82	88	6.9	0.21
Calprest	50	11	17	38	149	4	82	88	82	88	6.9	0.21
Phadia EliA	50	12	19	<15	183	12	75	95	90	86	14.3	0.26
	15						83	84	77	89	5.3	0.20
Certest	50	12	19	<50	>200	/	83	84	77	89	5.3	0.20

^aIf median FC was outside the measurement range, the detection limit was used to calculate the ratio IBD/non-IBD.

FC, faecal calprotectin; IBD, inflammatory bowel disease; LR+, positive likelihood ratio; LR-, negative likelihood ratio; NPV, negative predictive value; PPV, positive predictive value; QB BM, Quantum Blue Bühlmann.

patients, respectively. This resulted in a prevalence of 39% IBD. In the 19 other patients, IBD could be excluded. No discrepancies between endoscopy and histology were observed. In our patients without IBD, the following other organic diseases were diagnosed: nonsteroidal anti-inflammatory drug enteropathy ($n=1$), adenoma ($n=1$), and infectious colitis ($n=1$).

For all tested assays, a cut off of 50 $\mu\text{g/g}$ faeces is recommended by the manufacturer, except for Calfast (70 $\mu\text{g/g}$). We performed a ROC curve analysis (data not shown) to determine the cut off that provided the highest sum of sensitivity and specificity. This led to an optimal cut off for Bühlmann Quantum Blue of 75 $\mu\text{g/g}$ and for Phadia EliA of 15 $\mu\text{g/g}$. For the other assays, the optimal cut offs were identical to the one from the manufacturer. Sensitivity, specificity, predictive values, and likelihood ratios at the calculated cut offs are summarized in Table 4.

At the optimal cut offs, the assays were clinically (positive/negative) discordant in only two of the 31 patients (6.5%) and these discrepancies were only found in the non-IBD patients (false-positives). None of the patients with other organic disease ($n=3$) had an elevated calprotectin level with any of the assays.

In addition, we determined the discriminatory power of the tests by dividing the median calprotectin value in the IBD patients by the median value in the non-IBD patients. As shown in Table 4, this resulted for Bühlmann Quantum Blue, Phadia EliA, Calprolab ELISA, Calprest ELISA, and Calfast in the following ratios: 14, 12, 10, 4, and 3, respectively.

Follow up. Data from patients with UC and CD were analysed separately. We added the data from the newly diagnosed IBD patients (from the subpopulation with

suspected IBD) to those from the known IBD patients to obtain a larger sample size. As such, 21 patients with CD and 22 patients with UC were included in the analysis. Results are shown in Figure 1.

None of the manufacturers proposes cut offs for identifying relapse in patients with IBD. Analogous to the evaluation of calprotectin in the diagnostic setting, we calculated optimal cut offs for detecting endoscopic disease activity by ROC analysis (data not shown).

Crohn's disease. Median SES-CD score in the subgroup with CD was 3 (range 0–15, IQR 1–7). We defined three subgroups according to the corresponding endoscopic disease activity. Endoscopic remission was defined as a SES-CD score of 0, mild endoscopic disease activity as SES-CD 1–4, and moderate-to-severe endoscopic disease activity as SES-CD >4. These groups consisted of five, nine, and seven patients, respectively. Of the patients with active inflammation, 10/16 (63%) had ileal disease, 3/16 (19%) colonic disease, and 3/16 (19%) ileocolonic disease. None of our patients in endoscopic remission had any histological disease activity. For all patients with endoscopic disease activity (SES-CD ≥ 1), pathologists reported the presence of active inflammation.

The results of the calprotectin assays for the assessment of endoscopic disease activity in CD are listed in Table 5. Optimal cut offs for the quantitative assays ranged from 29 to 89 $\mu\text{g/g}$. Performance of calprotectin in detecting mild disease activity (SES-CD ≤ 4) was low for all methods, with a sensitivity between 22 and 44%. For patients with a higher disease activity level (SES-CD >4), sensitivity increased for the different assays to 83–100%. All the assays displayed a specificity of 80%,




Figure 1. Faecal calprotectin measured by different assays in patients grouped by Mayo endoscopic subscore for ulcerative colitis (a) and simple endoscopic score for Crohn's disease (b).

Solid lines indicate the median values of each group; dashed lines indicate the optimal cut offs determined by ROC curve analysis. QB BM, Bühlmann Quantum Blue.

except for Certest which only reached a specificity of 40% at a cut off of 50 $\mu\text{g/g}$. This assay showed a specificity of 100% at the higher cut off of 200 $\mu\text{g/g}$, but this resulted in substantially lower sensitivities.

Ulcerative colitis. As for CD, patients with UC were divided into three subgroups: endoscopic remission (Mayo endoscopic subscore 0, $n=7$), mild endoscopic disease activity (Mayo 1, $n=7$), and moderate-to-severe endoscopic disease activity (Mayo 2 or 3, $n=8$). No discrepancies were found between endoscopy and

histology regarding the absence/presence of disease activity.

Table 6 summarizes the performance characteristics of the different tests in UC. Optimal cut offs for the quantitative assays ranged from 31 to 129 $\mu\text{g/g}$. A Mayo score of 2–3 could be differentiated from remission with a sensitivity of 100% by all the assays, whereas the sensitivities to distinguish Mayo 1 from remission were between 71 and 100%. The high sensitivities of Certest (at a cut off of 50 $\mu\text{g/g}$) and Calfast were associated with lower specificities (57%

Table 5. Performance of the calprotectin assays for detecting disease activity in Crohn's disease

Assay	Cut off ($\mu\text{g/g}$)	SES-CD 1-4 ($n=9$) vs. SES-CD 0 ($n=5$)						SES-CD >4 ($n=7$, except Calprest/Calfast $n=6$) vs. SES-CD 0 ($n=5$)					
		Sensitivity	Specificity	NPV	PPV	LR+	LR-	Sensitivity	Specificity	NPV	PPV	LR+	LR-
Calprolab	89	22	80	36	67	1.1	1.0	86	80	80	86	4.3	0.18
QB BM	77.5	44	80	44	80	2.2	0.7	100	80	100	88	5	0
Phadia EliA	28.5	22	80	36	67	1.1	1.0	86	80	80	86	4.2	0.21
Calfast	75.5	33	80	40	75	1.7	0.8	83	80	80	83	3.3	0.2
Calprest	46.5	44	80	44	80	2.2	0.7	100	80	100	86	5	0
Certest	50	33	40	25	50	0.56	1.7	100	40	100	70	1.7	0
	200	22	100	42	100	/	0.8	71	100	71	100	/	0.3

IBD, inflammatory bowel disease; LR+, positive likelihood ratio; LR-, negative likelihood ratio; NPV, negative predictive value; PPV, positive predictive value; QB BM, Quantum Blue Bühlmann; SES-CD, simple endoscopic score for Crohn's disease.

Table 6. Performance of the calprotectin assays for detecting disease activity in ulcerative colitis

Assay	Cut off ($\mu\text{g/g}$)	Mayo 1 ($n=7$) vs. Mayo 0 ($n=7$, except Calprest/Calfast $n=6$)						Mayo 2-3 ($n=8$) vs. Mayo 0 ($n=7$, except Calprest/Calfast $n=6$)					
		Sensitivity	Specificity	NPV	PPV	LR+	LR-	Sensitivity	Specificity	NPV	PPV	LR+	LR-
Calprolab	86	71	86	75	83	5	0.33	100	86	100	89	6.0	0
QB BM	129	71	86	75	83	5	0.33	100	86	100	89	6.0	0
Phadia	31	86	86	86	86	6	0.17	100	86	100	89	6.0	0
Calfast	50	100	67	100	88	3	0	100	67	100	80	3.0	0
Calprest	56	86	83	83	86	5	0.17	100	83	100	89	6.0	0
Certest	50	100	57	100	70	2.3	0	100	57	100	73	2.3	0
	200	57	86	67	80	4.0	0.5	75	86	75	86	5	0.3

LR+, positive likelihood ratio; LR-, negative likelihood ratio; Mayo, Mayo endoscopic subscore; NPV, negative predictive value; PPV, positive predictive value; QB BM, Quantum Blue Bühlmann.

and 67%, respectively) in comparison to the other assays.

Discussion

The aim of this study was to determine the correlation between the most widely used calprotectin assays and to compare their clinical performance for diagnosis and follow up of IBD.

In the first part of our study, we evaluated the quantitative correlation between the different tests. Correlation coefficients were moderate to excellent. However, when performing regression analysis, significant differences (up to 5-fold) were observed, even though manufacturers propose identical cut offs. For that reason, analysing the same sample with different calprotectin assays will result in very different values for FC.

These results are in line with a recent study from the UK National External Quality Assessment Service that

revealed up to 3.8-fold differences between ELISAs from different manufacturers.⁸ Studies comparing the Bühlmann Quantum Blue rapid test with an ELISA kit found correlation coefficients between 0.56 and 0.94. Slopes of the regression curves varied from 1.0 to 2.24 and intercepts from -11 to -40.⁹⁻¹²

Some of the variability in our study might have been caused by the different extraction devices and procedures, as was shown in the study of Whitehead et al.⁸ However, we also observed major quantitative differences between assays using the same extraction device.

In the second part of our study, we assessed the diagnostic accuracy of FC for the diagnosis of IBD among patients with a suspect clinical history. The optimal cut offs determined by ROC curve analysis resulted in highly comparable performance characteristics for all assays. A recent meta-analysis of prospective studies in adults found sensitivity, specificity, positive predictive value, and negative predictive value of 93, 96, 91, and 97%, respectively.⁴ The majority of the pooled

results were obtained with Phical or Calprest ELISA. In our study, a slightly lower performance was found. False-positive FC results have been attributed to several causes (e.g. infections, malignancies, and several drugs).⁴ However, in our two patients having a false-positive result with all assays, further investigation revealed no explanation for their abdominal symptoms or elevated calprotectin. In three patients, another organic abdominal disease than IBD was diagnosed, but they all had a low FC.

Two newly diagnosed IBD patients had FC levels below the limit of detection, with 3/5 and 5/5 assays. These two false-negative results were from patients with mild, strictly ileal CD. It has been shown in previous studies that this type of disease can result in low calprotectin levels despite detectable endoscopic inflammation.^{13–15} A false-negative calprotectin can also be caused by very liquid stool samples.¹⁶ In this regard, we rejected one extremely aqueous stool sample in which the calprotectin concentration was below the detection limit for all assays, thereby affecting the performance of all assays equally. Although the influence of the bowel preparation on the first stool sample is expected to be minimal, a small dilution effect on the other samples cannot be completely excluded. This could have slightly lowered the sensitivity, negative predictive value, and optimal cut offs of all calprotectin tests to the same extent.

Besides specificity and sensitivity, we determined the calprotectin value in IBD patients relative to non-IBD patients for each of the investigated assays. The highest discriminatory power between IBD and non-IBD samples was achieved with Bühlmann Quantum Blue, Phadia EliA, and Calprolab ELISA, with ratios that were ≥ 10 .

In the last part of our study, we compared the different assays in their ability to detect endoscopic disease activity in the follow up of IBD patients. In CD, mild endoscopic disease activity (SES-CD ≤ 4) could not be reliably differentiated from remission with any of the assays. Sixty-six per cent (6/9) of these patients only showed mild ileal disease, which results in low calprotectin levels, as already discussed. However, a SES-CD score >4 could be identified with all the assays, except Certest, with good sensitivities. In UC, sensitivities for mild disease were moderate to good and substantially better than in mild Crohn's disease. All assays had a very good sensitivity in detecting moderate-to-severe UC.

Nine studies on the performance of calprotectin in the follow up of IBD patients were reviewed by Lewis et al.,² sensitivity and specificity ranged from 70 to 93% and from 70 to 92%, respectively. Another recent prospective study conducted in our centre, in which samples were sent to another laboratory for FC

measurement by PhiCal ELISA, included a large cohort of CD patients ($n=87$). Sensitivities and specificities of 60 and 79%, respectively, for a SES-CD score of ≥ 2 were observed. In the same study, 39 patients with UC were included, with sensitivities and specificities of 71 and 100% for Mayo scores of 1–3.¹⁷ Our performance characteristics are thus in line with these observations.

Our optimal cut offs for detecting endoscopic disease activity in IBD varied widely from one assay to another. In general, they were rather low compared to other publications. In the review by Lewis et al.,² cut offs between 50 and 200 $\mu\text{g/g}$ were used by the different authors. The study by D'Haens et al.¹⁷ proposed a cut off of 250 $\mu\text{g/g}$. We hypothesize that mean disease activity of our study population was lower than in other studies, and more patients had only very mild inflammation. However, our (sub)populations were rather small, which is the main limitation of our study. Therefore, more extensive studies are needed on the different FC assays, and their cut offs for follow up in IBD have to be established in larger cohorts.

Overall, the evaluated immunochromatography tests are easier and faster to perform than the ELISA tests and are acceptable alternatives to an ELISA for diagnostic determinations. For purposes of follow up, a quantitative test with a high measuring range is preferred. The Certest is only semiquantitative, whereas the two quantitative rapid tests in our study only have a low measurement range. However, Bühlmann recently developed a high range version of the Quantum Blue assay, with an extended measurement range (100–1800 $\mu\text{g/g}$), specifically for follow up samples.

Since the immunochromatographic fast tests have short turnaround times and can be performed on one individual sample at the time, they would be practical and cost-effective in laboratories with small amounts of samples. ELISA is more suited for large laboratories, where large amounts of samples can be analysed in batch, which is more time- and cost-effective.

The Phadia automated EliA is a new method that may combine the fast and easy manipulation of the rapid tests with the precision and high measurement range of an ELISA test. However, the measurement results and optimal cut offs that were obtained in this study were remarkably lower than those from other assays. The cut off for diagnosis was optimal at a level of 15 $\mu\text{g/g}$. This is as low as the detection limit of the assay, which is analytically unacceptable. Further studies with larger study populations are needed to evaluate the performance of this new assay and to establish its cut offs.

In conclusion, we found a acceptable performance of all calprotectin tests. However, large quantitative

differences between the assays were observed. This variation makes it impossible to use methods interchangeably and enforces the urgent need for a further standardization.

Acknowledgements

Calprotectin test kits were provided for free or at reduced prices by all manufacturers. The authors would like to thank Bühlmann, Phadia/ThermoFisher, Eurospital, Calpro, and Biotest for this support.

Funding

This research received no specific grant from any funding agency in the public, commercial, or not-for-profit sectors.

Conflict of interest

The authors declare that there is no conflict of interest.

References

1. Bernstein CN, Fried M, Krabshuis JH, et al. World Gastroenterology Organization practice guidelines for the diagnosis and management of IBD in 2010. *Inflamm Bowel Dis* 2010; 16: 112–114.
2. Lewis JD. The utility of biomarkers in the diagnosis and therapy of inflammatory bowel disease. *Gastroenterology* 2011; 140: 1817–1826.
3. Centre for Evidence-based Purchasing. Evidence review: value of calprotectin in screening out irritable bowel syndrome, CEP09026. London: NHS Purchasing and Supply Agency, 2009.
4. Van Rheenen PF, Van de Vijver E and Fidler V. Faecal calprotectin for screening of patients with suspected inflammatory bowel disease: diagnostic meta-analysis. *BMJ* 2010; 341: c3369.
5. De Vos M, Dewit O, D'Haens G, et al. Fast and sharp decrease in calprotectin predicts remission by infliximab in anti-TNF naïve patients with ulcerative colitis. *J Crohns Colitis* 2012; 6: 557–562.
6. Schroeder KW, Tremaine WJ and Ilstrup DM. Coated oral 5-aminosalicylic acid therapy for mildly to moderately active ulcerative colitis. A randomized study. *N Engl J Med* 1987; 317: 1625–1629.
7. Daperno M, D'Haens G, Van Assche G, et al. Development and validation of a new, simplified endoscopic activity score for Crohn's disease: the SES-CD. *Gastrointest Endosc* 2004; 60: 505–512.
8. Whitehead S, French J, Brookes M, et al. Between-assay variability of faecal calprotectin enzyme-linked immunosorbent assay kits. *Ann Clin Biochem* 2013; 50: 53–61.
9. Kolho KL, Turner D, Veereman-Wauters G, et al. Rapid test for fecal calprotectin levels in children with Crohn's disease. *J Pediatr Gastroenterol Nutr* 2012; 55: 436–439.
10. Coorevits L, Baert F and Vanpoucke H. Faecal calprotectin: comparative study of the Quantum Blue rapid test and an established ELISA method. *Clin Chem Lab Med* 2012; 51: 825–831.
11. Dolci A and Panteghini M. Comparative study of a new quantitative rapid test with an established ELISA method for faecal calprotectin. *Clin Chim Acta* 2012; 413: 350–351.
12. Wassell J, Wallage M and Brewer E. Evaluation of the Quantum Blue rapid test for faecal calprotectin. *Ann Clin Biochem* 2012; 49: 55–58.
13. Sipponen T, Savilahti E, Kolho KL, et al. Crohn's disease activity assessed by fecal calprotectin and lactoferrin: correlation with Crohn's disease activity index and endoscopic findings. *Inflamm Bowel Dis* 2008; 14: 40–46.
14. Sipponen T, Karkkainen P, Savilahti E, et al. Correlation of faecal calprotectin and lactoferrin with an endoscopic score for Crohn's disease and histological findings. *Aliment Pharmacol Ther* 2008; 28: 1221–1229.
15. Jones J, Loftus EV, Pannaccione R, et al. Relationships between disease activity and serum and fecal biomarkers in patients with Crohn's disease. *Clin Gastroenterol Hepatol* 2008; 6: 1218–1224.
16. Vestergaard TA, Nielsen SL, Dahlerup JF, et al. Fecal calprotectin: assessment of a rapid test. *Scand J Clin Lab Invest* 2008; 68: 343–347.
17. D'Haens G, Ferrante M, Vermeire S, et al. Fecal calprotectin is a surrogate marker for endoscopic lesions in inflammatory bowel disease. *Inflamm Bowel Dis* 2012; 18: 2218–2224.